

SECOND HAND BED CHECKLIST

If you've sold or bought/obtained an Icare Homecare Bed second hand, Icare Medical Group takes no responsibility of the product or transaction involved. Warranties are non-transferable. This form is designed for user peace of mind, that the bed has been sold/purchased in good faith that it is in good working order. This pre-purchase Second Hand Bed Checklist ensures all the key electrics and parts are in good working, and may prevent faults before the transaction is made. It is recommended to have a technician perform a maintenance inspection before the bed is sold, and after the bed has been installed in new home.

SELLERS RESPONSIBILITIES:

1: Before bed is advertised for sale.
Is the bed still covered by warranty?
If Yes; have the bed serviced by the dealer who provided the equipment. Faulty parts will be supplied free of charge to the dealer.
If No; have the bed serviced by the dealer who provided the equipment. Faulty parts will have to be purchased.

2. Record the Serial Number.
Located at the head of bed on a cross-bar near the control box. e.g. BPNS-A00000

3. Check over the bed.

Upholstery clean and free from tears	
All cables are in tact, no cuts, no exposed wires, firmly plugged in	
Handpiece lights are working in full	
Wipe all metal parts with disinfectant	
Completed maintenance tasks in User Manual and Handbook	

4. Seller understands the recommendations and Icare Medical Group does not take any responsibility of this transaction.

Signature

BUYERS RESPONSIBILITIES:

1: Find out about the bed.

How old is the bed	
Has the bed had regular maintenance and inspections?	
Has the bed had a recent service?	
Is the bed clean, upholstery in good condition, free from tears/fading	
Has the seller completed the 'Second Hand Bed Checklist'?	

2. Installation & Set-up.
Follow the instructions outlined in bed Handbook and User Manual.

3. Follow the Routine Maintenance Points.
Complete maintenance tasks outlined in the User Manual and Handbook.

4. Record the details of your nearest dealer.
In the event of your second hand bed having a fault, it is recommended you have contact details of a store near you to contact. Fees may apply for their services as warranties are non-transferable.

5. Buyer understands the recommendations and Icare Medical Group does not take any responsibility of this transaction.

Signature